

4 INCH

1.2 INCH

1 INCH

1.5 INCH

1.5 INCH

5 INCH

5 INCH

**APPLICATION INSTRUCTIONS:**

Gently buff the nail surface to remove keratin, rinse, and keep nails clean and dry. Choose the right nail size and apply glue or adhesive tape inside the selected nail. Attach the nail from the base to the tip, pressing firmly for 8 seconds.

**REMOVAL INSTRUCTIONS:**

**For Glue:** Soak nails in remover or wrap with aluminum foil for 10 minutes. Once softened, gently lift from the base and sides using a cuticle pusher.  
**For Adhesive Tape:** No remover needed. Gently lift the nail from both sides with a cuticle pusher.

**PRECAUTIONS:**

Avoid applying on inflamed or red areas. Keep out of children's reach. Do not force removal to prevent nail damage.

# PRESS ON NAILS

24 PCS OF BEAUTIFUL NAILS

4 INCH

0.5 INCH